



**Model KL CO<sub>2</sub> Fractional Laser**

- ▶ 30W glass tube



**Model KL-R CO<sub>2</sub> Fractional Laser**

- ▶ 30W/60W RF tube from top international suppliers



**Model KL CO<sub>2</sub> Laser**

- ▶ Single/repeat/continuous/ultra pulse
- ▶ Longer arm is easy for more operations



**CO<sub>2</sub> Laser Microsurgery System**

- ▶ CO<sub>2</sub> laser with micromanipulator for E.N.T. treatment



**Model JC-100D CO<sub>2</sub> Laser**

- ▶ Single/repeat/continuous/ultra pulse



**Model ZYS Helium-Neon Laser**

- ▶ It is for anti-inflammation, analgesia and promoting tissue growth



**CONTACT US**




**Model KL-M(H) Q Switched Nd:YAG Laser**

- ▶ 1064nm+532nm/Q switched mode and 300μs long pulse mode
- ▶ High energy/spot size is adjustable/flat cap spot



**Model KL-M(S) 900 Picosecond Laser**

- ▶ 900ps~2.5ns/4ns/4+4ns
- ▶ 1064nm+532nm



**Model KL-M(P) Picosecond Laser**

- ▶ 350ps/1.43GW high peak power



**Model KL-S**

**Series of LED Therapy Instrument**

- ▶ Optional red/yellow/blue light source



**Handheld LED Therapy Instrument**

- ▶ Red light source/rechargeable



**Model KL-L Series of IPL Therapy Instrument**

- ▶ Optional filters wavelength:420nm,515nm, 560nm, 590nm, 615nm, 640nm, 695nm, ACNE,VASCULAR,Calibration
- ▶ Professional mode/stacking mode SOPT technology



**Model KL-D Nd:YAG & Alexandrite Laser**

- ▶ 1064nm+755nm long-pulse laser with adjustable pulsewidth and spot size



**Model KL Nd:YAG Laser**

- ▶ 1064nm surgery laser
- ▶ Max.output power:60W
- ▶ Optional φ0.4mm or φ0.6mm fiber



**Model KL-P Pulsed Nd:YAG Laser**

- ▶ 1064nm surgery laser
- ▶ Max. energy 900mJ
- ▶ Optional φ0.4mm or φ0.6mm fiber



**CONTACT US**



*Model KL-600/KL-600C Diode Laser*

- ▶ 808nm hair removal laser
- ▶ Professional/sliding/stacking mode



*Model KL-300 Diode Laser*

- ▶ 808nm fiber coupled hair removal laser
- ▶ Long life span/light weight handpiece



*Model KL-E Erbium Fractional Laser*

- ▶ 1550nm fiber laser
- ▶ Fractional/sliding handpiece



*Model KL-Q1 Lithotripsy Laser*

- ▶ 1064nm Nd:YAG laser
- ▶ Effective and safe for treating kidney stone and gallstone



*Model GLK Erbium Fractional Laser*

- ▶ 2940nm/high quality RF tube
- ▶ Accurate/safe/low thermal damage



*Model KY-15D40 Air Cooling Machine*

- ▶ Variable air speed at 9 levels
- ▶ Continuous working or working at regular time

### Parts



Articulated arm



Laser module



Power supply



**CONTACT US**

Multimagotrade.com  
+40-758-667-671